

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinders, 250 ml, A, conf. cert.

Batch: 22-17




Serial:

Measuring Nr. 16497

Batch size:

1

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Measured value:	49,9268 ml	150,0765 ml	250,3626 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.12.2017

Balance:

Sartorius - LA
1200 S

DE-M