

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinders, 250 ml, A, conf. cert.




Batch: 14-18

Serial:

Measuring Nr. 19817

Batch size: 500

Volume:	50,000 ml	150,000 ml	250,000 ml
Tolerance:	± 1,000 ml	± 1,000 ml	± 1,000 ml
Average value:	49,8894 ml	150,0183 ml	250,0214 ml
Deviation:	0,2290 ml	0,1919 ml	0,1592 ml
Precision (CV):	0,4580 %	0,1279 %	0,0637 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	11.04.2018
Balance:	Sartorius - LA 1200 S

