

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinders, 250:2 ml, A, conf. cert.




Batch: 21-19

Serial:

Measuring Nr. 30211

Batch size: 1200

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	49,8219 ml	149,9147 ml	250,0041 ml
Deviation:	0,1393 ml	0,1806 ml	0,1773 ml
Precision (CV):	0,2787 %	0,1204 %	0,0709 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	29.05.2019
Balance:	Sartorius - LE244S-OCE

DE-M