

Certificate of performance



Article: 3 715 025

Graduated cylinder, 250:2 ml, A, conf. cert.




Batch: 19-21

Serial:

Measuring Nr.: 46796

Batch size: 500

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	50,0455 ml	150,0810 ml	250,0867 ml
Deviation:	0,2265 ml	0,1651 ml	0,1482 ml
Precision (CV):	0,4529 %	0,1101 %	0,0593 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	17.05.2021
Balance:	Sartorius - LA 1200 S

