

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 025

**Graduated cylinder, 250:2 ml, A, conf. cert.**




**Batch:** 34-24

**Serial:**

**Measuring Nr.:** 68868

**Batch size:** 500

Volume:	50,000 ml	150,000 ml	250,000 ml
<b>Tolerance:</b>	± 1,000 ml	± 1,000 ml	± 1,000 ml
<b>Average value:</b>	50,0985 ml	150,2179 ml	250,2472 ml
<b>Deviation:</b>	0,2403 ml	0,1761 ml	0,1635 ml
<b>Precision (CV):</b>	0,4807 %	0,1174 %	0,0654 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	27.08.2024
<b>Balance:</b>	Sartorius - LA 1200 S

