

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 025

Graduated cylinder, 250:2 ml, A, conf. cert.




Batch: 49-25

Serial:

Measuring Nr. 75761

Batch size: 699

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
Tolerance:	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
Average value:	50,0901 ml	150,1156 ml	250,1681 ml
Deviation:	0,2388 ml	0,1793 ml	0,1898 ml
Precision (CV):	0,4776 %	0,1195 %	0,0759 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	09.12.2025
Balance:	Sartorius - LA 1200 S

