

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 BL

Grad. cylinder, 500:5 ml, blue grad., A, conf. cert.




Batch: 36-23

Serial:

Measuring Nr. 64409

Batch size: 500

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	99,9333 ml	300,0081 ml	500,3247 ml
Deviation:	0,2058 ml	0,1326 ml	0,2074 ml
Precision (CV):	0,2058 %	0,0442 %	0,0415 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	13.09.2023
Balance:	Sartorius - LA 1200 S

