

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 050 BL

**Grad. cylinder, 500:5 ml, blue grad., A, conf. cert.**




**Batch:** 33-25

**Serial:**

**Measuring Nr.** 74743

**Batch size:** 398

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Average value:	100,2483 ml	300,4641 ml	500,2990 ml
Deviation:	0,2595 ml	0,2814 ml	0,2236 ml
Precision (CV):	0,2595 %	0,0938 %	0,0447 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

29.09.2025

**Balance:**

Sartorius -  
MSE14202S-  
OCE-D0

