

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050 USP

Grad. cylinder, 500:5 ml, USP, A, conf. cert.




Batch: 33-25

Serial:

Measuring Nr. 74174

Batch size: 400

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,0000 ml	± 2,0000 ml	± 2,0000 ml
Average value:	99,8198 ml	300,3783 ml	500,5107 ml
Deviation:	0,4710 ml	0,4178 ml	0,3647 ml
Precision (CV):	0,4710 %	0,1393 %	0,0729 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	14.08.2025
Balance:	Sartorius - LA 1200 S

