

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 050

Graduated cylinders, 500 ml, A, conf. cert.

Batch: 33-17




Serial:

Measuring Nr. 16498

Batch size:

1

Volume:	100,000 ml	300,000 ml	500,000 ml
Tolerance:	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
Measured value:	100,4203 ml	300,1432 ml	499,8712 ml



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.12.2017

Balance:

Sartorius - LA
1200 S

DE-M