

# Certificate of performance



**Article:** 3 715 050

**Graduated cylinders, 500:5 ml, A, conf. cert.**




**Batch:** 20-18

**Serial:**

**Measuring Nr.:** 20794

**Batch size:** 500

Volume:	100,000 ml	300,000 ml	500,000 ml
<b>Tolerance:</b>	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
<b>Average value:</b>	100,2194 ml	300,2382 ml	500,1824 ml
<b>Deviation:</b>	0,5421 ml	0,4140 ml	0,3954 ml
<b>Precision (CV):</b>	0,5421 %	0,1380 %	0,0791 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	44
<b>Date:</b>	18.05.2018
<b>Balance:</b>	Sartorius - LA 1200 S

