

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 BL

Grad. cylinder, 1000:10 ml, blue grad., A, conf. cert.




Batch: 24-26

Serial:

Measuring Nr. 78665

Batch size: 100

Volume:	200,0000 ml	600,0000 ml	1000,0000 ml
Tolerance:	$\pm 5,0000$ ml	$\pm 5,0000$ ml	$\pm 5,0000$ ml
Average value:	200,4248 ml	600,6704 ml	1000,3941 ml
Deviation:	0,5333 ml	0,8199 ml	0,5005 ml
Precision (CV):	0,2666 %	0,1367 %	0,0501 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

22.06.2026

Balance:

Sartorius -
MSE14202S-
OCE-D0

DE-M