

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 USP

Graduated cylinder, 1000 ml, USP, A, conf. cert.




Batch: 31-24

Serial:

Measuring Nr.: 68380

Batch size: 90

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 3,0000 ml	± 3,0000 ml	± 3,0000 ml
Average value:	200,6354 ml	599,9338 ml	999,4540 ml
Deviation:	0,5184 ml	0,4620 ml	0,5004 ml
Precision (CV):	0,2592 %	0,0770 %	0,0500 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	02.08.2024
Balance:	Sartorius - MSE14202S- OCE-D0

