

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100 USP

Graduated cylinder, 1000 ml, USP, A, conf. cert.




Batch: 24-26

Serial:

Measuring Nr. 78644

Batch size: 150

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	$\pm 3,0000$ ml	$\pm 3,0000$ ml	$\pm 3,0000$ ml
Average value:	199,4492 ml	600,0807 ml	1000,9513 ml
Deviation:	0,6912 ml	0,8808 ml	0,6557 ml
Precision (CV):	0,3456 %	0,1468 %	0,0656 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

18.06.2026

Balance:

Sartorius -
MSE14202S-
OCE-D0

DE-M