

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100

Graduated cylinder, 1000:10 ml, A, conf. cert.




Batch: 11-22

Serial:

Measuring Nr. 53316

Batch size: 599

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	199,8116 ml	599,6119 ml	999,4271 ml
Deviation:	0,7678 ml	0,6033 ml	0,9422 ml
Precision (CV):	0,3839 %	0,1006 %	0,0942 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

23.03.2022

Balance:

Sartorius - LC
6200 S

DE-M