

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 100

Graduated cylinder, 1000:10 ml, A, conf. cert.




Batch: 26-23

Serial:

Measuring Nr. 63589

Batch size: 500

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,4942 ml	599,8411 ml	999,6062 ml
Deviation:	0,7596 ml	0,6477 ml	0,7530 ml
Precision (CV):	0,3798 %	0,1080 %	0,0753 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

30.06.2023

Balance:

Sartorius -
MSE14202S-
OCE-D0

