

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 100

**Graduated cylinder, 1000:10 ml, A, conf. cert.**




**Batch:** 03-26

**Serial:**

**Measuring Nr.:** 76192

**Batch size:** 368

Volume:	200,000 ml	600,000 ml	1000,000 ml
Tolerance:	± 5,0000 ml	± 5,0000 ml	± 5,0000 ml
Average value:	200,2909 ml	599,7990 ml	998,6887 ml
Deviation:	0,9599 ml	0,6003 ml	0,5698 ml
Precision (CV):	0,4799 %	0,1000 %	0,0570 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	4
Date:	20.01.2026
Balance:	Sartorius - MSE14202S- OCE-D0

