

Certificate of performance



Article: 3 715 200 BL

Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.




Batch: 31-18

Serial:

Measuring Nr.: 23171

Batch size: 41

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	402,2714 ml	1201,6618 ml	1999,2978 ml
Deviation:	1,4646 ml	1,5370 ml	2,2124 ml
Precision (CV):	0,3662 %	0,1281 %	0,1106 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	44
Date:	08.08.2018
Balance:	Sartorius - LC 6200 S

