

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 200 BL

**Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.**




**Batch:** 26-19

**Serial:**

**Measuring Nr.:** 30897

**Batch size:** 50

Volume:	400,000 ml	1200,000 ml	2000,000 ml
<b>Tolerance:</b>	± 10,000 ml	± 10,000 ml	± 10,000 ml
<b>Average value:</b>	401,8198 ml	1199,1739 ml	1997,4354 ml
<b>Deviation:</b>	1,4028 ml	1,3571 ml	0,9780 ml
<b>Precision (CV):</b>	0,3507 %	0,1131 %	0,0489 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	25.06.2019
<b>Balance:</b>	Sartorius - LC 6200 S

