

Certificate of performance



Article: 3 715 200 BL

Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.




Batch: 30-22

Serial:

Measuring Nr.: 57675

Batch size: 53

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,0000 ml	± 10,0000 ml	± 10,0000 ml
Average value:	397,4966 ml	1196,8764 ml	1998,8747 ml
Deviation:	1,4083 ml	1,1624 ml	0,9515 ml
Precision (CV):	0,3521 %	0,0969 %	0,0476 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	27.07.2022
Balance:	Sartorius - LC 6200 S

