

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 200 BL

Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.

Batch: 32-24




Serial:

Measuring Nr.: 68537

Batch size:

250

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	397,9807 ml	1198,2120 ml	1999,4857 ml
Deviation:	1,0332 ml	1,0335 ml	1,4616 ml
Precision (CV):	0,2583 %	0,0861 %	0,0731 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	13.08.2024
Balance:	Sartorius - MSE14202S- OCE-D0

