

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 200 BL

**Grad. cylinder, 2000:20 ml, blue grad., A, conf. cert.**




**Batch:** 16-26

**Serial:**

**Measuring Nr.** 77531

**Batch size:** 38

Volume:	400,000 ml	1200,000 ml	2000,000 ml
<b>Tolerance:</b>	± 10,000 ml	± 10,000 ml	± 10,000 ml
<b>Average value:</b>	401,3918 ml	1200,8991 ml	1999,8071 ml
<b>Deviation:</b>	1,3112 ml	1,5271 ml	2,0672 ml
<b>Precision (CV):</b>	0,3278 %	0,1273 %	0,1034 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	16.04.2026
<b>Balance:</b>	Sartorius - MSE14202S- OCE-D0

