

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 200 USP

Graduated cylinder, 2000:20 ml, USP, A, conf. cert.




Batch: 36-17

Serial:

Measuring Nr. 14072

Batch size: 24

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 6,0000 ml	± 6,0000 ml	± 6,0000 ml
Average value:	400,2439 ml	1199,7860 ml	1998,7287 ml
Deviation:	2,2056 ml	0,4565 ml	0,3947 ml
Precision (CV):	0,5514 %	0,0380 %	0,0197 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	12.09.2017
Balance:	Sartorius - LC 6200 S

