

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 200

**Graduated cylinder, 2000:20 ml, A, conf. cert.**




**Batch:** 26-19

**Serial:**

**Measuring Nr.** 31179

**Batch size:** 150

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	400,0141 ml	1199,9771 ml	2000,6057 ml
Deviation:	0,2409 ml	0,4322 ml	0,6163 ml
Precision (CV):	0,0602 %	0,0360 %	0,0308 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:	7
Date:	02.07.2019
Balance:	Sartorius - LE244S-OCE

