

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 715 200

Graduated cylinder, 2000:20 ml, A, conf. cert.




Batch: 47-20

Serial:

Measuring Nr.: 42762

Batch size: 87

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,0000 ml	± 10,0000 ml	± 10,0000 ml
Average value:	399,7354 ml	1199,8686 ml	2000,5139 ml
Deviation:	0,2947 ml	0,6134 ml	0,6144 ml
Precision (CV):	0,0737 %	0,0511 %	0,0307 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.11.2020
Balance:	Sartorius - LC 6200 S

DE-M