

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 715 200

**Graduated cylinder, 2000:20 ml, A, conf. cert.**




**Batch:** 45-25

**Serial:**

**Measuring Nr.:** 75201

**Batch size:** 140

Volume:	400,000 ml	1200,000 ml	2000,000 ml
<b>Tolerance:</b>	± 10,000 ml	± 10,000 ml	± 10,000 ml
<b>Average value:</b>	398,2995 ml	1199,0643 ml	2000,1065 ml
<b>Deviation:</b>	2,0467 ml	2,2867 ml	2,2286 ml
<b>Precision (CV):</b>	0,5117 %	0,1906 %	0,1114 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	06.11.2025
<b>Balance:</b>	Sartorius - MSE14202S- OCE-D0

