

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 001

Grad. cylinders, 10 ml, NS 10, A, conf. cert.




Batch: 03-18

Serial:

Measuring Nr. 17174

Batch size: 125

Volume:	2,0000 ml	6,0000 ml	10,0000 ml
Tolerance:	± 0,1000 ml	± 0,1000 ml	± 0,1000 ml
Average value:	1,9819 ml	5,9829 ml	9,9837 ml
Deviation:	0,0141 ml	0,0084 ml	0,0100 ml
Precision (CV):	0,7045 %	0,1404 %	0,1004 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	9
Date:	18.01.2018
Balance:	Sartorius - LE244S-OCE

