

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 005

Grad. cylinder, 50:1 ml, NS 19, A, conf. cert.




Batch: 42-24

Serial:

Measuring Nr.: 69933

Batch size: 200

Volume:	10,000 ml	30,000 ml	50,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	9,9934 ml	29,9971 ml	49,9825 ml
Deviation:	0,0520 ml	0,0368 ml	0,0518 ml
Precision (CV):	0,5203 %	0,1226 %	0,1036 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	04.11.2024
Balance:	Sartorius - MSE224S-0CE -DU

