

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 010 BL

Grad. cylinders, 100:1 ml, NS 24, blue, A, conf. cert.

Batch: 06-26




Serial:

Measuring Nr. 76834

Batch size:

241

Volume:	20,000 ml	60,000 ml	100,000 ml
Tolerance:	± 0,5000 ml	± 0,5000 ml	± 0,5000 ml
Average value:	20,0553 ml	60,0842 ml	100,0742 ml
Deviation:	0,1078 ml	0,0903 ml	0,0786 ml
Precision (CV):	0,5391 %	0,1506 %	0,0786 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	27.02.2026
Balance:	Sartorius - LA 1200 S

