

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 025 BL

**Grad. cylinder, 250:2 ml, NS 29, blue grad., A, conf. cert**




**Batch:** 10-26

**Serial:**

**Measuring Nr.:** 77157

**Batch size:** 50

Volume:	50,0000 ml	150,0000 ml	250,0000 ml
<b>Tolerance:</b>	± 1,0000 ml	± 1,0000 ml	± 1,0000 ml
<b>Average value:</b>	50,2536 ml	150,1357 ml	250,1273 ml
<b>Deviation:</b>	0,0973 ml	0,0628 ml	0,1444 ml
<b>Precision (CV):</b>	0,1946 %	0,0419 %	0,0578 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	4
<b>Date:</b>	26.03.2026
<b>Balance:</b>	Sartorius - LA 1200 S

