

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 754 050 BL

**Grad. cylinder, 500:5 ml, blue, NS 34, A, conf. cert.**




**Batch:** 30-19

**Serial:**

**Measuring Nr.** 31997

**Batch size:** 90

Volume:	100,000 ml	300,000 ml	500,000 ml
<b>Tolerance:</b>	± 2,5000 ml	± 2,5000 ml	± 2,5000 ml
<b>Average value:</b>	100,1127 ml	300,1260 ml	500,1930 ml
<b>Deviation:</b>	0,4667 ml	0,2715 ml	0,6514 ml
<b>Precision (CV):</b>	0,4667 %	0,0905 %	0,1303 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

<b>Operator:</b>	9
<b>Date:</b>	06.08.2019
<b>Balance:</b>	Sartorius - LE244S-OCE

