

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 754 200

Grad. cylinder, 2000 ml, NS 45, A, conf. cert.




Batch: 48-17

Serial:

Measuring Nr. 16528

Batch size: 25

Volume:	400,000 ml	1200,000 ml	2000,000 ml
Tolerance:	± 10,000 ml	± 10,000 ml	± 10,000 ml
Average value:	400,9012 ml	1200,2110 ml	1999,0156 ml
Deviation:	2,3649 ml	1,6703 ml	1,1061 ml
Precision (CV):	0,5912 %	0,1392 %	0,0553 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 4788

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	16.12.2017
Balance:	Sartorius - LC 6200 S

