

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 00114

Volumetric flasks, 10 ml, NS 14, A, conf. cert.

Batch: 17-26

Serial:

Measuring Nr. 77764

Batch size: 447


Volume: 10,0000 ml

Tolerance: $\pm 0,0600$ ml

Average value: 10,0055 ml

Deviation: 0,0116 ml

Precision (CV): 0,1164 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

27.04.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M