

# Certificate of performance



ISO 9001  
CERTIFIED

**Article:** 3 662 002

**Volumetric flasks, 20 ml, NS 10, A, conf. cert.**

**Batch:** 17-26

**Serial:**

**Measuring Nr.** 78195

**Batch size:** 1177


**Volume:** 20,0000 ml

**Tolerance:**  $\pm 0,0400$  ml

**Average value:** 20,0048 ml

**Deviation:** 0,0047 ml

**Precision (CV):** 0,0233 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

7

**Date:**

01.06.2026

**Balance:**

Sartorius -  
LE244S-OCE

DE-M