

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 0030 14

Volumetric flask, 30 ml, NS 14, A, conf. cert.,

Batch: 50-21


Serial:

Measuring Nr. 51104

Batch size:

100

Volume:	30,000 ml
Tolerance:	$\pm 0,0600$ ml
Average value:	30,0103 ml
Deviation:	0,0067 ml
Precision (CV):	0,0222 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

14.12.2021

Balance:

Sartorius -
LE244S-OCE

DE-M