

# Certificate of performance

**Article:** 3 662 005 USP

**Volumetric flask, 50 ml, NS 12, USP, A, conf. cert.**


**Batch:** 13-26

**Serial:**

**Measuring Nr.** 77169

**Batch size:** 489

<b>Volume:</b>	<b>50,000 ml</b>
<b>Tolerance:</b>	$\pm 0,0500$ ml
<b>Average value:</b>	49,9993 ml
<b>Deviation:</b>	0,0134 ml
<b>Precision (CV):</b>	0,0267 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

**witeg Labortechnik GmbH**

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

**Operator:**

4

**Date:**

30.03.2026

**Balance:**

Sartorius -  
MSE224S-0CE  
-DU