

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 006 FB

Vol. flask, 50 ml, NS 14, flat bottom, A, conf. cert


Batch: 09-22

Serial:

Measuring Nr. 53508

Batch size: 150

Volume:	50,0000 ml
Tolerance:	$\pm 0,1000$ ml
Average value:	50,0032 ml
Deviation:	0,0053 ml
Precision (CV):	0,0106 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	24.03.2022
Balance:	Sartorius - LE244S-OCE

