

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 011

Volumetric flask, 100 ml, NS 14, A, conf. cert.

Batch: 16-26

Serial:

Measuring Nr. 77747

Batch size: 1189


Volume: 100,0000 ml

Tolerance: $\pm 0,1000$ ml

Average value: 100,0123 ml

Deviation: 0,0166 ml

Precision (CV): 0,0166 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

24.04.2026

Balance:

Sartorius -
MSE224S-0CE
-DU

DE-M