

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 025 BL USP

Vol. flasks, 250 ml, NS 14, blue, USP, A, conf. cert.


Batch: 27-21

Serial:

Measuring Nr. 47944

Batch size: 277

| | |
|------------------------|-----------------|
| Volume: | 250,000 ml |
| Tolerance: | $\pm 0,1200$ ml |
| Average value: | 249,9970 ml |
| Deviation: | 0,0163 ml |
| Precision (CV): | 0,0065 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

| | |
|------------------|--------------------------|
| Operator: | 4 |
| Date: | 16.07.2021 |
| Balance: | Sartorius - LA 1200 S |

DE-M