

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 050 W

Vol. flask., 500 ml, NS 19, white grad., A, conf.cert.,

Batch: 51-20


Serial:

Measuring Nr. 43509

Batch size:

145

Volume:	500,000 ml
Tolerance:	$\pm 0,2500$ ml
Average value:	500,0611 ml
Deviation:	0,0330 ml
Precision (CV):	0,0066 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	12.01.2021
Balance:	Sartorius - LC 6200 S

DE-M