

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 060 29 BL

Volumetric flask, 600 ml, NS 29, blue, class A, conf. cer


Batch: 44-24

Serial:

Measuring Nr. 69874

Batch size: 10

Volume:	600,000 ml
Tolerance:	$\pm 0,6000$ ml
Average value:	600,1513 ml
Deviation:	0,0264 ml
Precision (CV):	0,0044 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	4
Date:	30.10.2024
Balance:	Sartorius - LA 1200 S

