

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 125

Volumetric flasks, 125 ml, NS 14, A, conf. cert.


Batch: 17-23

Serial:

Measuring Nr. 62767

Batch size: 50

Volume:	125,000 ml
Tolerance:	$\pm 0,1000$ ml
Average value:	125,0117 ml
Deviation:	0,0192 ml
Precision (CV):	0,0154 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator: 7

Date: 26.04.2023

Balance: Sartorius -
LE244S-OCE

DE-M