

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 662 300

Volumetric flasks, 3000 ml, NS 29, A, conf. cert.

Batch: 05-19


Serial:

Measuring Nr. 26948

Batch size:

32

Volume:	3000,0000 ml
Tolerance:	$\pm 0,8000$ ml
Average value:	3000,3781 ml
Deviation:	0,1332 ml
Precision (CV):	0,0044 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

22.01.2019

Balance:

Sartorius -
MSE14202S-
OCE-DO

DE-M