

# Certificate of performance



ISO 9001  
CERTIFIED

Article: 3 662 400 BL

Vol. flask, 4000 ml, NS 34, blue grad., A, conf. cert.

Batch: 15-26


Serial:

Measuring Nr. 77588

Batch size:

15

|                 |                 |
|-----------------|-----------------|
| Volume:         | 4000,0000 ml    |
| Tolerance:      | $\pm 1,2000$ ml |
| Average value:  | 4000,8108 ml    |
| Deviation:      | 0,1484 ml       |
| Precision (CV): | 0,0037 %        |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16  
97877 Wertheim  
Germany  
+49 (0) 9342 / 9301-0  
info@witeg.de

Operator:

4

Date:

21.04.2026

Balance:

Sartorius -  
MSE14202S-  
OCE-D0

DE-M