

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 000

Volumetric flask, 1 ml, NS 7, trapezoidal, A, conf. cert.

Batch: 46-18


Serial:

Measuring Nr. 25634

Batch size:

592

| | |
|------------------------|-----------------|
| Volume: | 1,0000 ml |
| Tolerance: | $\pm 0,0250$ ml |
| Average value: | 1,0023 ml |
| Deviation: | 0,0041 ml |
| Precision (CV): | 0,4124 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

4

Date:

14.11.2018

Balance:

Sartorius -
LE244S-OCE

DE-M