

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 669 002 BL

Vol. flasks, 3 ml, NS 7, trap., blue grad., A, conf. cert.

Batch: 30-23


Serial:

Measuring Nr. 63874

Batch size:

149

Volume:	3,0000 ml
Tolerance:	$\pm 0,0250$ ml
Average value:	3,0008 ml
Deviation:	0,0045 ml
Precision (CV):	0,1492 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

28.07.2023

Balance:

Sartorius -
LE244S-OCE

DE-M