

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 002 USP

Vol.flask, 20 ml, USP, NS 10, amber glass, A, conf.cert


Batch: 41-25

Serial:

Measuring Nr. 75005

Batch size: 497

Volume:	20,000 ml
Tolerance:	$\pm 0,0300$ ml
Average value:	19,9967 ml
Deviation:	0,0027 ml
Precision (CV):	0,0134 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:	7
Date:	22.10.2025
Balance:	Sartorius - LE244S-OCE

DE-M