

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 002 USP

Vol.flask, 20 ml, USP, NS 10, amber glass, A, conf.cert

Batch: 19-26

Serial:

Measuring Nr. 78082

Batch size:

267


Volume: 20,0000 ml

Tolerance: $\pm 0,0300$ ml

Average value: 20,0015 ml

Deviation: 0,0037 ml

Precision (CV): 0,0187 %



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

7

Date:

21.05.2026

Balance:

Sartorius -
LE244S-OCE

DE-M