


Certificate of performance

Article number 3 670 003
Batch 05-17
Serial number
Measuringnr./Batch size 10082 / 1200
Balance Sartorius - LE244S-0CE
Operator 4
Date 09.02.2017

Volumetric flask 25 ml NS 10 Class A KB brown glass

Volume 25,0000ml
Tolerance 0,0400ml
Average volume 24,9933 ml
Deviation 0,0060 ml
Number of measurings 35



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
| info@witeg.de

