

Certificate of performance



ISO 9001
CERTIFIED

Article: 3 670 003

Vol. flasks, 25 ml, NS 10, amber glass, A, conf. cert.

Batch: 24-20


Serial:

Measuring Nr. 38871

Batch size:

1192

| | |
|------------------------|--------------------|
| Volume: | 25,0000 ml |
| Tolerance: | ± 0,0400 ml |
| Average value: | 24,9962 ml |
| Deviation: | 0,0095 ml |
| Precision (CV): | 0,0379 % |



These final test values refer to the delivered nominal volume (max. volume), instrument and liquid (dist. water) at equilibrium with ambient temperature (20 °C). The test was effected according to ISO 1042

witeg Labortechnik GmbH

Am Bildacker 16
97877 Wertheim
Germany
+49 (0) 9342 / 9301-0
info@witeg.de

Operator:

9

Date:

30.06.2020

Balance:

Sartorius -
LE244S-OCE

DE-M